ReDura™

Instructions for Use Regenerative Dural Repair Patch

Contents

Instructions for Use Regenerative Dural Repair Patch (ReDura™).....	P1
Gebrauchsanweisung Regenerative Dural Repair Patch(ReDura™)	P5
Instrucciones de uso Parche regenerativo de reparación dural (ReDura™)	P10
(ReDura™) لصقة إصلاح الجافية التجديدية تعليمات الاستخدام	P14

Instructions for Use Regenerative Dural Repair Patch


Description:


Regenerative Dural Repair Patch owns a 3D micro-fibrous network resembling the micro-structure of human dural matrices which provides a temporary scaffold for migration and growth of dural cells and tissues. As the material degrades gradually in human body, the new dural tissue regenerates and the defect is repaired. In three months after product implantation, the newborn tissue can form relatively stable repair substitute. The intended life time of the product is 2 years, considering the product is absorbed in body within 2 years.

Performance Characteristics:

Regenerative Dural Repair Patch is an absorbable synthetic patch with biomimetic design. This patch offers good biocompatibility, high strength and excellent handling as well as superior cerebrospinal fluid (CSF) leak resistance. The device can be used onlay or sutured. It can be used in combination with water soluble glues, such as: fibrin glue, PEG sealant, etc.

Models & Specifications: RDS-X

Shape	Model	Dimensions(±5mm)		Units/Case	Thickness (mm)
Square/ rectangle	RDS-1	15×20		1	0.1~0.5
	RDS-2	20×30		1	
	RDS-3	30×40		1	
	RDS-4	40×60		1	
	RDS-5	60×60		1	
	RDS-6	60×80		1	
	RDS-7	60×140		1	
	RDS-8	80×80		1	
	RDS-9	80×120		1	
	RDS-10	100×150		1	
	RDS-11	150×150		1	
	RDS-12	25×25(1 inch × 1 inch)		1	
	RDS-13	25×75(1 inch × 3 inch)		1	
	RDS-14	50×50(2 inch × 2 inch)		1	
	RDS-15	75×75(3 inch × 3 inch)		1	
	RDS-16	100×125(4 inch × 5 inch)		1	
	RDS-1205	25×25(1 inch × 1 inch)		5	
	RDS-1305	25×75(1 inch × 3 inch)		5	
RDS-1405	50×50(2 inch × 2 inch)		5		
RDS-1505	75×75(3 inch × 3 inch)		5		
RDS-17	30×60		1		
RDS-18	30×90		1		
Shape (Oval)	Model	Specifications			
		a (±5mm)	b (±5mm)	Thickness (mm)	
	RDS-19	75	45	0.1~0.5	
	RDS-20	100	85		

Shape (Fan-shaped)	Model	Specifications				
		a (±5mm)	c (±5mm)	b (±5mm)	Central angle	Thickness (mm)
	RDS-21	105	35	67	60°	0.1~0.5

Intended Use:

This product is applied in surgically repairing operations in case of destructive dural defects.

Indications:

Regenerative Dural Repair Patch is used for repairing dural defects resulting from craniocerebral trauma, brain tumor, cerebrovascular disease, congenital disease of nervous system, posterior fossa surgery, spinal canal disease and other dura mater defects caused in various neurosurgical procedures and accidents.

Patient Target Groups:

The product can be used for patients with dura mater defects resulting from craniocerebral trauma, brain tumor, cerebrovascular disease, congenital disease of nervous system, posterior fossa surgery, spinal canal disease and other dura mater defects caused in various neurosurgical procedures and accidents.

The safety and effectiveness of Regenerative Dural Repair Patch has been studied in adult patients with no restriction of gender. The patients should not have the history of severe allergy and serious immunodeficiency and their health status need to meet the standard requirements for the neurosurgical procedures.

Intended User& Environment:

The product must be used by clinicians, and they must carefully read the instruction before use.

Operation in surgery room with professional health care facilities. Normative and local regulations must be observed.

Contraindications:

1. This product is forbidden in patients with open craniocerebral injury and low-immunity.
2. Not recommended to cover dural mater defects involving mastoid air cells.

Adverse Effects:

Possible complications that may occur with any dural repair material include fever, infection, adhesion, cerebrospinal fluid leaks, inflammation, irritation, neurologic deficit, meningitis, pseudomeningocele, brain edema, hemorrhage, hydrocephalus, abnormal thickening, graft failure, re-operation and discomfort in the body.

Expected Clinical Benefits:

The clinical benefits of ReDura™ include but not limit:

1. ReDura™ can degrade gradually in human body, then the new dural tissue regenerates.
2. ReDura™ can effectively repair the dural mater defects and reduce the incidence of CSF leakage.
3. ReDura™ is made of synthetic material, which has advantage in no risk of infecting the virus; easy handling; good conformability; good mechanical properties like strength and resistance to traction.

Material and Substances:

The material and substances can be exposed to patients, including:
Poly (L-lactide) (≤1.3g/piece);

Lactic acid (≤0.043g/d);

1,1,1,3,3,3-Hexafluoro-2-propanol (≤100 µg/g) ;

Ethyl alcohol (≤300 µg/g);

Shelf Life:

Shelf Life: 36 months.

Storage:

Stored at room temperature (5°C-35°C), kept from lights, and sealed.

Sterilization Method:

Sterilized by gamma radiation and supplied in sterilized condition.

Instructions for Surgical Implantation:

1. Inspect the package carefully before use.
2. Tear or cut the outer package along the notch and take out the inner package by aseptic manipulation.
3. The aseptic inner package can be placed directly in sterilized environment.
4. Cover the product on the dural defect and fix it by a routine suturing. During suturing, pinholes should stay on the positions 2-3mm away from the product patch edge. It is recommended that Regenerative Dural Repair Patch can be used together with water soluble glue, i.e., Fibrin glue, PEG sealant, etc.
5. Alternatively, Regenerative Dural Repair Patch can also be used as an onlay graft by trimming it to a size permitting a minimum of 15-20mm overlap with the remaining dura. Stay sutures may be used if desired. For best results, the Regenerative Dural Repair Patch should be used as described below.
Step 1: Place the article into an appropriately sized container of saline.
Step 2: Press onto the article with gloved fingers to ensure it absorbs the saline (2-5 minutes).
Step 3: Remove the article from the saline once it becomes transparent to aseptically cut into appropriately sized and shaped sections. Keep Regenerative Dural Repair Patch immersed in sterile saline until ready for use.
6. This product can be trimmed and cut as needed. The used product must be disposed as medical waste instead of being sterilized for reuse.
7. The user could close skull after being stitched up and make sure that there is no CSF leakage.
8. Mild folding does not affect the quality and performance of the product.

Warning and Precautions:

- This sterilized product shall be used immediately once unpacked.
- The clinician shall select a suitable-size device based on the shape and size of the patient defect. It is recommended that size of the patch is at least 15-20mm larger than that of the defect.
- Do not use the device if the package is damaged or there is any visible damage on the patch.
- The device is only for single use. Re-use of the device might pose a risk of infection to patients.
- If the device is to be sutured, any commercially available suture can be used, and tensionless suturing technique must be used to prevent tearing the device.
- If the glue is used together, the commercially available water-soluble glues can be used, such as: fibrin glue, PEG sealant, etc.
- Never use the device with non-water-soluble glues.
- Do not re-sterilize it.
- Do not engage in any strenuous physical activity in half a month after surgery.

Explanation on graphics, symbols, abbreviations for labels:

 : DO NOT RE-USE	 : DO NOT RSTERILIZE
 : CAUTION	STERILE R : STERILIZED USING IRRADIATION
LOT : BATCH CODE	 : CONSULT INSTRUCTIONS FOR USE OR CONSULT ELECTRONIC INSTRUCTIONS FOR USE
SN : SERIAL NUMBER	 : USE-BY DATE
REF : CATALOGUE NUMBER	 : CE MARKING OF CONFORMITY AND IDENTIFICATION NUMBER OF NOTIFIED BODY
 : DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE	 : TEMPERATURE LIMIT
 : KEEP DRY	 : KEEP AWAY FROM SUNLIGHT
 : FRAGILE, HANDLE WITH CARE	MR : MR SAFE
MD : MEDICAL DEVICE	UDI : UNIQUE DEVICE IDENTIFIER
 : SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE	 : DATE OF MANUFACTURE
EC REP : AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/ EUROPEAN UNION	 : IMPORTER
 : MANUFACTURER	

The Link of Summary of Safety and Clinical Performance: (It has not yet been generated until the SSCP is uploaded to EUDAMED).

For the electronic version of instruction for use (e-IFU), please refer to the website: www.medprin.com.

 **Medprin Regenerative Medical Technologies Co., Ltd.**
3/F & 7/F, Zone E, 80 Lanyue Road, Science City, New High-Tech
Industrial Park, Guangzhou 510663, China

Customers Service:
Sales and Marketing, Medprin Regenerative Medical Technologies Co., Ltd.
TEL: + 86-20-32296118 FAX: +86-20-32296128

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Any serious incident related to this device, please report to the manufacturer (Email: ae@medprin.com) and the competent authority of the Member State.

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 **Medprin Regenerative Medical Technologies Co., Ltd.**
3/F & 7/F, Zone E, 80 Lanyue Road, Science City, New High-Tech Industrial Park,
Guangzhou 510663, China

Kundendienst:
Vertrieb und Marketing, Medprin Regenerative Medical Technologies Co, Ltd.
TEL: + 86-20-32296118 FAX: +86-20-32296128

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