NeoDura™

Instructions for Use
Dural Repair Patch

Contents

Instructions for Use Dural Repair Patch (NeoDura™)	P1
لصقة الجافية التصحيحية تعليمات الاستخدام (NeoDura™)	P5
Gebrauchsanweisungen Dural Repair Patch (NeoDura™)	P9
Instrucciones de Uso Parche (matriz recuperadora) de reparación dural (NeoDura™).....	P13
Mode d'emploi Patch de Réparation Durale (NeoDura™)	P17
Istruzioni per L'Uso Patch Durale (NeoDura™)	P21
Instruções de Uso Adesivo de Reparação Dural (NeoDura™)	P25
Руководство по эксплуатации Пластырь для восстановления твёрдой мозговой оболочки (NeoDura™)	P29
Kullanım Talimatları Dural Onarım Yaması (NeoDura™)	P33

Instructions for Use

Dural Repair Patch

Description

Dural Repair Patch owns a 3D micro-fibrous network resembling the micro-structure of human dural matrices which provides a temporary scaffold for migration and growth of dural cells and tissues. As the material degrades gradually in human body, the new dural tissue regenerates and the defect is repaired. In three months after product implantation, the newborn tissue can form relatively stable repair substitute. The intended life time of the product is 1 year, considering the product is absorbed in body within 1 year.

Performance Characteristics

Dural Repair Patch is manufactured with absorbable material polylactic acid and porcine gelatin by emerging biomimetic technology, the materials are widely used in clinic for decades. This onlay patch offers extremely high strength and excellent handling and conforming capabilities as well as superior cerebrospinal fluid (CSF) leak resistance. The device has multi-use: onlay and suturable. It can be used in combination with water-soluble glues, such as: Fibrin glue, PEG sealant, etc.

Models and Specifications

Model	Dimensions (width × length) (±0.5cm)	Equivalent English Dimensions	Units/Case	Thickness (mm)
RDP-1	2.5 cm × 2.5 cm	(1 inch × 1 inch)	1	0.1-0.5
RDP-105	2.5 cm × 2.5 cm	(1 inch × 1 inch)	5	
RDP-2	2.5 cm × 7.5 cm	(1 inch × 3 inch)	1	
RDP-205	2.5 cm × 7.5 cm	(1 inch × 3 inch)	5	
RDP-3	3 cm × 4 cm	-	1	
RDP-4	4 cm × 6 cm	-	1	
RDP-5	5 cm × 5 cm	(2 inch × 2 inch)	1	
RDP-505	5 cm × 5 cm	(2 inch × 2 inch)	5	
RDP-6	6 cm × 6 cm	-	1	
RDP-7	6 cm × 8 cm	-	1	
RDP-8	7.5 cm × 7.5 cm	(3 inch × 3 inch)	1	
RDP-805	7.5 cm × 7.5 cm	(3 inch × 3 inch)	5	
RDP-9	6 cm × 14 cm	-	1	
RDP-10	8 cm × 12 cm	-	1	
RDP-11	10 cm × 12.5 cm	(4 inch × 5 inch)	1	
RDP-12	15 cm × 15 cm	(6 inch × 6 inch)	1	

Intended Use

This product is used as a dural repair patch for the surgical repair of dural defects.

Indications

Dural Repair Patch is used for repairing dural defects resulting from craniocerebral trauma, brain tumor, cerebrovascular disease, congenital disease of nervous system, posterior fossa surgery, spinal canal disease and other dura mater defects caused in various neurosurgical procedures and accidents.

Patient Target Groups

The product can be used for patients with dura mater defects resulting from craniocerebral trauma, brain tumor, cerebrovascular disease, congenital disease of nervous system, posterior fossa surgery, spinal canal disease and other dura mater defects caused in various neurosurgical procedures and accidents.

The safety and effectiveness of Dural Repair Patch has been studied in adult patients with no restriction of gender. The patients should not have the history of severe allergy and serious immunodeficiency and their health status need to meet the standard requirements for the neurosurgical procedures.

Intended user/ environment

The product must be used by clinicians, and they must carefully read the instruction before use.

Operation in surgery room with professional health care facilities. Normative and local regulations must be observed.

Contraindications

1. Dural Repair Patch is forbidden for use in any other conditions besides the indications described above.
2. This product contains highly purified derivative gelatin from porcine skin, so it should not be used in patients with a known history of hypersensitivity to porcine gelatin materials.
3. Dural Repair Patch should not be used in patients with unresolved serious infections.
4. Not recommended to cover dural defects involving mastoid air cells.
5. This product is forbidden in patients with open craniocerebral injury and low-immunity.

Warnings

1. This product is sterile medical device for single use, and must not be re-sterilized or reused. Otherwise it will increase the risk of infection and cerebrospinal fluid leakage.
2. Do not use the product if the original package is damaged. Any package damage will destruct the sterile barrier system, causing contamination to the product or perhaps failure of its function. The product should be used immediately once opened.
3. No clinical experience is available of using Dural Repair Patch in pregnant women.

Precautions

1. The clinician can select the suitable model based on the size of the defects in patients. It is recommended that the size of the patch is at least 15-20mm larger than that of the defect.
2. Rinse surgical gloves to remove any glove powder prior to handling the product.
3. Only physicians qualified in the appropriate surgical techniques can use this product.
4. It is the responsibility of the physician to inform the patient of the porcine gelatin content of the patch before use.
5. Dural Repair Patch should be cut as needed. The remaining material after cutting should be disposed as medical waste.

6. If the device is to be sutured, any commercially available suture can be used and tensionless suturing technique must be used to prevent tearing the device.
7. If the glue is used together, the commercially available water-soluble glues can be used, such as: Fibrin glue, PEG sealant, etc. Don't use it with the non-water-soluble glues.
8. Do not engage in any strenuous physical activity in half a month after surgery.
9. Occurrence of serious non-microbial tissue response may require removal of the implant.
10. If an infection develops, it should be treated aggressively. Unresolved infection may require removal of the patch.

Adverse Effects

Possible complications that may occur with any dural repair material include fever, infection, adhesion and cerebrospinal fluid leaks, inflammation, irritation, neurologic deficit, meningitis, pseudomeningocele, brain edema, hemorrhage, hydrocephalus, abnormal thickening, graft failure, re-operation and discomfort in the body.

Expected Clinical Benefits

The clinical benefits of NeoDura™ include but not limit:

1. NeoDura™ can degrade gradually in human body, then the new dural tissue regenerates.
2. NeoDura™ can effectively repair the dural mater defects and reduce the incidence of CSF leakage.
3. NeoDura™ is made of polylactic acid and gelatin, it has good biological performance with tissue; easy handling, good conformability; good mechanical properties like strength, and resistance to traction.

Material and Substances

The material and substances can be exposed to patients, including:

- Poly-L-lactic acid (≤1.04g/piece)
- Gelatin (≤0.325g/piece)
- Lactic acid (≤34.7 mg/d)
- 1,1,1,3,3,3-Hexafluoro-2-propanol (≤100µg/g)
- Ethyl alcohol (≤300µg/g)

Directions for Use

1. Inspect the package and read the Instructions for Use carefully prior to use.
2. Tear the aluminum-plastic package and take out the primary package by sterile operation. The primary package is sterile and can be placed directly in a sterilized environment.
3. Cut dry Dural Repair Patch to fit the size of the defect, ensuring the patch can cover over 15-20mm of the patient's remaining dura mater around the edge.
4. Soak the patch in adequate amount of normal saline, and 3-5 minutes (with mild shaking) is recommended for best conformity. After hemostasis, place the soaked patch directly on the defect by aseptic technique. Dural Repair Patch can also be applied in the dry state where necessary, by directly placing on the defect by aseptic technique and moistening with saline. It is recommended that Dural Repair Patch can be used together with water-soluble glue, i.e., Fibrin glue, PEG sealant, etc.

5. Dural Repair Patch could also be anchored by suturing attributing to its high strength. This suturable feature is beneficial to the repair of large defects and other conditions where suturing is required. During suturing, pinholes should stay on the positions 2-3mm away from the product patch edge. Absorbable suture is recommended for suturing.
6. CSF leakage must be observed after implantation and it should be confirmed that there is no leakage prior to skull closure.

Sterilization

Dural Repair Patch is sterilized by gamma radiation and supplied as sterile.








Shelf Life

36 months.

Storage

This product should be stored in clean and ventilated environment between 5°C and 30°C, free from corrosive gas.

Symbols used on labelling

	: Do not re-use		: Do not re-sterilize
	: Caution		: Sterilized using irradiation
	: Batch code		: Consult instructions for use or consult electronic instruction for use
	: Serial number		: Use-by date
	: Catalogue number		: Do not use if package is damaged and consult instruction for use
	: Fragile, handle with care		: The temperature limit is 5°C to 30°C.
	: Manufacturer		: CE-marking of conformity and identification number of Notified Body
	: MR safe		: Medical Device
	: Unique Device Identifier		: Contains biological material of animal origin
	: Single sterile barrier system with protective packaging inside		: Date of manufacture

The Link of Summary of Safety and Clinical Performance (SSCP)

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI (4260333310002PG). The URL to the Eudamed public website is <https://ec.europa.eu/tools/eudamed>.

For the electronic version of instruction for use (e-IFU), please refer to the website: www.medprin.com.

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Any serious incident related to this device, please report to the manufacturer (Email: ae@medprin.com) and the competent authority of the Member State.

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**Medprin Biotech GmbH**
Satış ve Pazarlama, Medprin Biotech GmbH